

DATA INPUT KEY

The following data points are required for the ClinicalQ Report Writer

LEGEND

CT: Central regions (Cz)

EO: Eyes Open

EC: Eyes Closed

UT: Under Task (e.g., reading)

TB: Theta/Beta ratio

LOCATION Cz

CT1: % change EO to EC Alpha

CT2: % change EO (before EC) to EO (after EC)

CT3: EO Theta/Beta

CT4: UT Theta/Beta

CT5: % change T/B EO to T/B UT

CT6: Total Amplitude

CT7: Theta/SMR

LOCATION O1

OC1: % change in Alpha EO to EC

OC2: % change EO (before EC) to EO (after EC)

OC3: Theta/Beta EO

OC4: Theta/Beta EC

OC5: % change T/B EO to T/B EC

LOCATIONS F3 and F4

FF1: % absolute difference between F3 and F4 Beta

FF2: % absolute difference between F3 and F4 Alpha

FF3: % absolute difference between F3 and F4 Theta

FF4: F3 Theta/Alpha amplitude ratio

FF5: F4 Theta/Alpha amplitude ratio

FF6: F3 Theta/Beta amplitude ratio

FF7: F4 Theta/Beta amplitude ratio

FF8: % absolute difference between F3 and F4 T/B

LOCATION Fz

FZ1: Delta amplitude

FZ2: HiBeta&Gamma/Beta

FZ3: Sum of HiBeta&Gamma & Beta

FZ4: LoAlpha/HiAlpha amplitude ratio

FREQUENCY RANGES USED FOR THE CLINICAL Q ASSESSMENT

Delta 1.5–2.5 Hz

Theta 3–7 Hz

Low Alpha 8–9 Hz

High Alpha 11–12 Hz

Alpha 8–12 Hz

SMR 13–15 Hz

Beta 16–25 Hz

High Beta-Gamma 28–40 Hz